

NEAR VISION AREA WITHOUT BINOCULAR DISTORTIONS IN PERSONALIZED FREE-FORM PROGRESSIVE LENSES WITH AN INCREASING CURVATURE FRONTAL BASE CURVE

Gago C*, Chamorro E*, Pascual E*, Subero M*, Cleva JM*, Alonso J*, Crespo D*
*Indizen Optical Technologies S.L., Spain

Corresponding author Melisa Subero

iot


e200

INTRODUCTION

The ideal optical **base curve** for each dioptic is represented in the **Tscherning Ellipse**: a graphical representation of the frontal surface power according to the total lens power in the best-form lenses (Figure 1).

The **Tscherning ellipse** shows two possible solutions:

- **The Oswalt branch**: represents the flatter curve most appropriate for a lens, which has become the basis for standardized lens base curve.
- **The Wollaston branch**: the most curved frontal surface.



When the power of a lens is selected with its **ideal base curve**, the **wearer enjoys a clear vision with minimal oblique astigmatism**. On the other hand, if the base curve falls outside the ideal range for a determined optical power, the wearer's off-axis visual acuity decreases precipitously. This performance has a great influence in a progressive lens since the **reading zone can only be accessed in an oblique angle**. Therefore, if the base curve is poorly matched, the patient might experience an abrupt drop-off in visual acuity of the reading area.

A **progressive addition lens (PAL)** has many optical powers; and therefore has many ideal base curves, which increase a diopter from the top to the bottom. The **distance visual area needs a flatter base curve**, whereas a **steeper base curve is most appropriate for the near visual area**. When a free-form progressive lens is manufactured with a single vision lens blank, various optical powers must share a single base curve. Consequently, the near zone does not have a correspondent frontal surface. In order to improve the base selection on the single vision lens blank, a continuously increasing base curve (three diopters from the top of the lens blank to the bottom) is proposed to have a better suited progressive prescription.

PURPOSE

To evaluate the differences in near vision areas without binocular distortion in personalized free-form progressive lenses due to the use of an increasing variable frontal base curve of the lens.

NEAR VISION AREA WITHOUT BINOCULAR DISTORTIONS IN PERSONALIZED FREE-FORM PROGRESSIVE LENSES WITH AN INCREASING CURVATURE FRONTAL BASE CURVE

Gago C*, Chamorro E*, Pascual E*, Subero M*, Cleva JM*, Alonso J*, Crespo D*

*Indizen Optical Technologies S.L., Spain

Corresponding author Melisa Subero

iot

e200

MATERIALS AND METHODS

DESIGN



A double-blind prospective observational longitudinal trial was performed on a sample of 21 presbyopic subjects (15 women and 6 men between 50 -63 years old), who have been worn progressive lenses for at least one year.

CLINICAL TRIAL PROCEDURE

After wearing a pair of tested PALs for a week, the undistorted near visual area was evaluated using a special test from IOT.

The test consists of a chart with different lines of letters of visual acuity 0.1 LogMAR, which patients have to choose their favorite lecture line (Figure 2) at 35 cm, without moving their heads.

Afterwards, the selected line is separated from the rest (Figure 3), and patients read the letters from the central one to their right/left areas until the blurry vision starts (**clear near vision area**) and/or they cannot recognize any more characters (**limit near vision area**).



STATISTICAL ANALYSIS

Design of randomized complete block test was used to determine differences in visual acuity among the lens designs. All statistical tests were performed using Statgraphics Centurion XVI.II software.

NEAR VISION AREA WITHOUT BINOCULAR DISTORTIONS IN PERSONALIZED FREE-FORM PROGRESSIVE LENSES WITH AN INCREASING CURVATURE FRONTAL BASE CURVE

Gago C*, Chamorro E*, Pascual E*, Subero M*, Cleva JM*, Alonso J*, Crespo D*

*Indizen Optical Technologies S.L., Spain


Corresponding author Melisa Subero




MATERIALS AND METHODS *continued*

PROGRESSIVE LENSES

2 progressive addition lenses, were tested for a week each. Both lenses, which were developed by IOT, showed the same power distribution perceived by the users, according to Sheedy's research (Figure 4).



However, one of the PALs showed an **aspheric base curve**, which increased its curvature from the top to the bottom part of the lens (Camber™ lens blank, Younger Optics). The other PAL had a simple **spherical base curve** (Figure 5).



NEAR VISION AREA WITHOUT BINOCULAR DISTORTIONS IN PERSONALIZED FREE-FORM PROGRESSIVE LENSES WITH AN INCREASING CURVATURE FRONTAL BASE CURVE

Gago C*, Chamorro E*, Pascual E*, Subero M*, Cleva JM*, Alonso J*, Crespo D*
*Indizen Optical Technologies S.L., Spain

Corresponding author Melisa Subero

iot


e200

RESULTS

The clinical assessment of the undistorted near vision area revealed significant differences among the different PALs tested. A **significantly wider and clearer near visual area was found in the progressive lenses with an increasing variable frontal base curve.**

A significant **increase of 9% in clear vision** of the near visual area was found in patients wearing lenses with an increasing variable frontal base curve (Figure 6).

UNDISTORTED NEAR VISION AREA



	PAL-A	PAL-B	P-VALUE
CLEAR NEAR VISION AREA (MM)	13.13 ± 3.28	14.32 ± 3.60	0.006*
LIMIT NEAR VISION AREA (MM)	17.17 ± 3.71	18.25 ± 3.57	0.070

Figure 6

CONCLUSIONS

The area of near vision without distortion can be improved using a variable frontal base curve that increases in curvature from top to bottom. Specifically, in this work, an increase of 9% in the clear area of near vision was found, indicating a better visual performance for reading tasks.

REFERENCES

- Han SC, Graham AD, Lin MC. (2011) Clinical Assessment of a Customized Free-Form Progressive Add Lens Spectacle. *Optom Vis Sci*, 88, 234-243.
- E. Sheedy J, Campbell C, King-Smith E, R. Hayes J. (2005) Progressive Powered Lenses: the Minkwitz Theorem. *Optometry and Vision Science*, 82(10), 1-9.
- Selenow A, A. Bauer E, R Ali S, Wayne Spencer L, J. Ciuffreda K. (2002) Assessing Visual Performance with Progressive Addition Lenses. *Optometry and Vision Science*, 79 (8), 502-505.
- Malacara Z, Malacara D. (1985) Tscherning ellipses and ray tracings in ophthalmic lenses. *American Journal of Optometry & Physiological Optics*, 62(7), 447 -455.
- Sheedy J, F. Raymond H., R. Hayes J. (2006) Progressive Addition Lenses – Measurements and ratings. *Optometry*, 77, 23-39.